

# Measurement of Leaf Gas Exchange Based on a New Single-Step CO<sub>2</sub> Response Method for Rapidly Obtaining A vs. C<sub>i</sub> Curves

John Ertle<sup>1</sup>, James Bunce<sup>2</sup>, Laura E. Dougherty<sup>2</sup>, Tim Doyle<sup>1</sup>, Graham Walter<sup>1</sup>  
<sup>1</sup> PP Systems, Amesbury, MA, <sup>2</sup> USDA-ARS, Beltsville, MD



## Introduction

- Ramping of CO<sub>2</sub> either up or down, and recording data once stabilized is a common technique used to develop A vs. C<sub>i</sub> curves.
- Researchers collecting A vs. C<sub>i</sub> curves may be limited by additional steps to validate and collect accurate data, including:
  - › Time-intensive (empty) chamber calibration prior to data collection.
  - › Post-processing of data after collection.
- We tested a novel, single-step CO<sub>2</sub> response (SSCO<sub>2</sub>R™) to eliminate extra calibration and data processing steps. The validation method was performed on leaves of tomato plants in two growing conditions.

## Materials & Methods

- A vs. C<sub>i</sub> data was captured for tomato leaves under field and growth chamber conditions using traditional steady-state methods and the single-step CO<sub>2</sub> response method (described below).
  - Both methods were tested on the same leaves for two tomato cultivars in the field, and three tomato cultivars in growth chambers.
  - All data was collected using the Ciras-4 Portable Photosynthesis System.



### The SSCO<sub>2</sub>R™ Method:

For the SSCO<sub>2</sub>R™ method, an adjustable volume plunger was inserted into the reference air line to reduce measurement differences in CO<sub>2</sub> between reference and analysis infrared gas exchange analyzers (IRGA).




### Data collection:

- A leaf cuvette was enclosed on a fully expanded tomato leaf.
- Steady state data collection was conducted for eight CO<sub>2</sub> concentrations between 100 and 700 μmol mol<sup>-1</sup>. Data was post-processed and an empty chamber ramp was conducted before collection on the leaf.
- SSCO<sub>2</sub>R™ data collection was conducted every one second during CO<sub>2</sub> ramping from 100 to 700 μmol mol<sup>-1</sup> at a rate of 200 μmol mol<sup>-1</sup> min<sup>-1</sup>.

## Results


Using the new single-step CO<sub>2</sub> response method, complete A vs. C<sub>i</sub> curves can be obtained for leaves of C<sub>3</sub> species in about 5 minutes, with no post-processing of data required.




**Figure 1.** Relationships between CO<sub>2</sub> assimilation rate (A) and internal CO<sub>2</sub> (C<sub>i</sub>) for a tomato leaf measured at 22 °C and 1500 mmol m<sup>-2</sup> s<sup>-1</sup> PPFD near midday in the field—either under steady-state CO<sub>2</sub> conditions or during an upward ramping of CO<sub>2</sub>—at a rate of 200 mmol mol<sup>-1</sup> min<sup>-1</sup>, using the SSCO<sub>2</sub>R™ Method.

## Results (continued)

- Traditional methods using steady-state determination of A vs C<sub>i</sub> was adequate to produce a curve for all tested cultivars in field and growth chamber conditions with eight total data points per curve. Collection of data (not including data processing and device setup) was roughly 30 minutes for each curve.
- The SSCO<sub>2</sub>R™ method was able to produce A vs C<sub>i</sub> curves with a high resolution of 300 data points per curve. Collection of data (not including device setup) was roughly five minutes for each curve.



**Figure 2.** CO<sub>2</sub> assimilation rate (A) and internal CO<sub>2</sub> (C<sub>i</sub>) for leaves of two different cultivars of tomato measured at 22 °C, and 1500 mmol m<sup>-2</sup> s<sup>-1</sup> PPFD within a few minutes of each other, near midday, in the field. Curves were determined using the SSCO<sub>2</sub>R™ Method with CO<sub>2</sub> increasing at a rate of 200 μmol mol<sup>-1</sup> min<sup>-1</sup>.



**Figure 3.** CO<sub>2</sub> assimilation rate (A) and internal CO<sub>2</sub> (C<sub>i</sub>) for leaves of three different lines of tomatoes grown in a growth chamber at 22 °C, a PPFD of 400 mmol m<sup>-2</sup> s<sup>-1</sup> for 12 hours per day, and measured at 22 °C, and 1500 mmol m<sup>-2</sup> s<sup>-1</sup> PPFD. Curves were determined using the SSCO<sub>2</sub>R™ Method with CO<sub>2</sub> increasing at a rate of 200 μmol mol<sup>-1</sup> min<sup>-1</sup>.

### Photosynthesis model parameters for tomato curves

Variable	Steady-state	SSCO <sub>2</sub> R™	Field Cultivar 1	Field Cultivar 2	Indoor Cultivar 1	Indoor Cultivar 2	Indoor Cultivar 3
V <sub>cmax</sub>	124	128	98	108	81	71	56
J <sub>max</sub>	205	209	189	220	174	161	120
TPU	10.1	9.9	9	12	10	9.5	6.1
Rd	3.8	4.3	3.3	2.9	2.3	3	3.6
g <sub>m</sub>	3.1	3.2	7.5	7	7	8.2	8

**Table 1.** Parameters of the Farquhar, von Caemmerer, Berry photosynthesis model derived from the A vs. C<sub>i</sub> curves of tomato leaves (Fig. 1-3) in this study. Units are mmol m<sup>-2</sup> s<sup>-1</sup> for V<sub>cmax</sub>, J<sub>max</sub>, TPU, and Rd, and μmol m<sup>-2</sup> s<sup>-1</sup> Pa<sup>-1</sup> for g<sub>m</sub>. For all measurements, leaf temperature was 22 °C, and PPFD was 1500 mmol m<sup>-2</sup> s<sup>-1</sup>.

## Conclusion

- Using the new single-step CO<sub>2</sub> response method, complete A vs. C<sub>i</sub> curves can be obtained for leaves of C<sub>3</sub> species in approximately five minutes with no post-processing of data required.
- Traditional methods for determining A vs C<sub>i</sub> curves took six times longer than the SSCO<sub>2</sub>R™ method.
- The SSCO<sub>2</sub>R™ method can reduce data collection time, affording researchers more time to take measurements without the need for data validation steps after collection.